



Development of near-field laser ablation inductively coupled plasma mass spectroscopy for sub-micrometric analysis of solid samples

C. Jabbour, J.-L. Lacour, M. Tabarant, A. Semerok, S. V. Fomichev, F. Chartier

► To cite this version:

C. Jabbour, J.-L. Lacour, M. Tabarant, A. Semerok, S. V. Fomichev, et al.. Development of near-field laser ablation inductively coupled plasma mass spectroscopy for sub-micrometric analysis of solid samples. 2016 Winter Conference on Plasma Spectrochemistry, Jan 2016, Tucson, United States.
hal-02445724

HAL Id: hal-02445724

<https://cea.hal.science/hal-02445724>

Submitted on 20 Jan 2020

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.

DE LA RECHERCHE À L'INDUSTRIE

**Development of near-field laser ablation
inductively coupled plasma mass
spectrometry for sub-micrometric
analysis of solid samples.**

C. JABBOUR¹, J-L. LACOUR¹, M. TABARANT¹,
A. SEMEROK¹, S.V. FOMICHEV², F. CHARTIER¹

1 CEA Saclay, DEN-DANS, Department of Physical Chemistry, Analytical Studies and Surface Reactivity Service, Nuclear Isotopic and Elemental Analytical development Laboratory, F-91191 Gif-sur-Yvette, France.

2 Kurchatov Institute, Moscow, Russia.
chirelle.jabbour@cea.fr

**2016 Winter Conference on Plasma Spectrochemistry,
January 10 - 16, 2016, Tucson, Arizona.**


SUMMARY

- ➡ *Objective*
- ➡ *Near Field Laser Ablation*
- ➡ *Experimental*
- ➡ *Multiparametric study: Results and Modeling Code*
- ➡ *Conclusions and Prospects*

OBJECTIVE




- Powerful analytical qualitative and quantitative method at ambient pressure.
- Optical diffraction limit → Lateral resolution μm




- Break the resolution limit → Sub-micrometer scale:
- High-resolution surface analysis.


NEAR FIELD LASER ABLATION




No tip enhancement effect
(No ablation)



Laser energy enhancement
(Near field ablation)



EXPERIMENTAL




Laser Nd:YAG (266 nm) , pulse duration 4 ns.

AFM: Atomic Force Microscope.

ICP-MS: Inductively Coupled Plasma - Mass Spectrometry : double focusing sector field

MULTIPARAMETRIC STUDY



Tip nature

Laser irradiance

Sample Nature

Tip dimensions

Laser wavelength

Distance tip-sample

Laser pulse duration

Environment gas



Parameters influencing the near field ablation efficiency, resolution, craters dimensions (amount of ablated mass)

RESULTS

Tip effect on the measurement resolution

$\lambda_{\text{laser}} = 266 \text{ nm}$; $\tau = 4 \text{ ns}$

- ❖ *Different tip natures and dimensions (diameter 10 - 30 nm → NO VISIBLE CRATERS).*
- ❖ *Observed craters only with the conductive diamond coated tip (Si + Diamond coating doped with B (diameter ~ 200 - 250 nm).*



The diameter of the obtained craters ranged from 100 to 120 nm using a single laser pulse.

RESULTS

Tip-sample distance and laser irradiance effects on the craters dimensions

$\lambda_{\text{laser}} = 266 \text{ nm}$; $\tau = 4 \text{ ns}$


- Tip-sample distance 

 - Crater diameter ~ 
 - Crater Depth 
 - Ablated mass 

- Laser Irradiance 

 - Crater diameter ~ 
 - Crater Depth 
 - Ablated mass 

Au (Single laser pulse)



RESULTS

Sample nature and multiple laser pulses effects on the craters dimensions


$\lambda_{\text{laser}} = 266 \text{ nm}$; $\tau = 4 \text{ ns}$

➤ Sample Nature

Samples: Gold / Silicon


- Crater diameter ~ 
- Crater depth : Au > Si
- ➡ Ablated mass: Au > Si

Sample: Gold
Irradiance: 80 mJ/cm^2
Tip-sample distance: 5 nm




➡ Ablated mass 

Number of laser pulses	Ablated mass	
	Au	Si
1	2 E+06 Atoms	6 E+05 Atoms
10	~ 4 fg	~ 0.3 fg
40	~ 12 fg	~ 1 fg
200	2 E+08 Atoms	4 E+07 Atoms




MODELING OF THE SURFACE HEATING

Si, 266 nm, $\alpha^{-1} = 5 \text{ nm}$, $\tau = 4 \text{ ns}$, $2r_0 = 1 \mu\text{m}$



Si, 266 nm, $\alpha^{-1} = 5 \text{ nm}$, $\tau = 4 \text{ ns}$, $2r_0 = 20 \text{ nm}$



τ : laser pulse duration

α^{-1} : sample absorption coefficient

$2r_0$: interaction size at the surface

$T = f(r)$ pour $z = 0$ et $t = \tau_{\text{laser}}$

Si (266nm): $\alpha^{-1} = 5 \text{ nm}$

$D_t = 0,8 \text{ cm}^2/\text{s}$

$L_t = 1,2 \mu\text{m}$

Au (266nm): $\alpha^{-1} = 11 \text{ nm}$


$D_t = 1,3 \text{ cm}^2/\text{s}$

$L_t = 1,4 \mu\text{m}$

MODELING OF THE SURFACE HEATING

Laser ablation on a large surface of interaction

Si, 266 nm, $\alpha^{-1} = 5 \text{ nm}$, $\tau = 4 \text{ ns}$, $2r_0 > 1 \mu\text{m}$




Ablation dominated by 1D sample thermal diffusion related to the pulse duration

Si (266nm): $\alpha^{-1} = 5 \text{ nm}$
 $D_t = 0,8 \text{ cm}^2/\text{s}$
 $L_t = 1,2 \mu\text{m}$

Laser ablation on a very small area of interaction (near field)

Si, 266 nm, $\alpha^{-1} = 5 \text{ nm}$, $\tau = 4 \text{ ns}$, $2r_0 = 20 \text{ nm}$




Ablation dominated by 3D sample thermal diffusion and absorption coefficient


Au (266nm): $\alpha^{-1} = 11 \text{ nm}$
 $D_t = 1,3 \text{ cm}^2/\text{s}$
 $L_t = 1,4 \mu\text{m}$

FIRST TESTS WITH ICP-MS

Ablated gold particles were collected on a filter using a simple aspiration tube set at a few mm from the tip before being analyzed by ICP-MS (25 000 craters /single laser pulse)



CONCLUSIONS AND PROSPECTS



- Near field laser ablation on gold sample and silicone wafer.
- Completion the multi-parametric study (wavelength, laser pulse duration...).
- Development of the AFM-ICP-MS interface (on-going).
- Particle characterization: particle size, mass, shape ... (on-going).
- Development of the elemental analysis method by ICP-MS.

Thank
you for
your
attention



MAY THE ATOMIC FORCE BE WITH ME !!

Commissariat à l'énergie atomique et aux énergies alternatives

Centre de Saclay | 91191 Gif-sur-Yvette Cedex

T. +33 (0)1 69 08 26 29 |

Etablissement public à caractère industriel et commercial | R.C.S Paris B 775
685 019

DEN-DANS

DPC

SEARS